
CEREBRO Head Assembly Parts Description Author: Dr. Vikram Pal Singh

Dept. of Psychology, UCSD Date: 01/06/23

Overall description:

- The head assembly is held in to the animal's head using two 3D printed horn that are permanently glued to the animal's headpost.
- 2. All the parts are printed with Titanium (Ti-6-Al-4V) using Direct Metal Laser Sintering (DLS).
- 3. All the 3D parts are held together with a 2 mm (diameter) stainless steel rod curved to match the profile of the animal's forehead (shown in yellow).



Detailed description of individual parts:

** I am using the same names here that you will find in the folder with 3D models.

1. New headpost.stl



Description:

Two of these headposts are glued on the monkey's headpost (using acrylic) while fitting the CEREBRO head assembly.

Material:

PP-Like Translucent White (Somos 9120)

2. Head Post cover.stl



Description:

These cover the headposts so that they do not break in the cage or get stuck in the mesh of the cage.

** It has a through hole on the side. I tap M4 screw though that and use a M4 set screw to hold this piece firmly on the headpost.

Material:

PP-Like Translucent White (Somos 9120)

CEREBRO Head Assembly Parts Description Author: Dr. Vikram Pal Singh

Dept. of Psychology, UCSD Date: 01/06/23

3. Head_Post_Holder.step	
	Description: There is a through hole along the short axis of this cylindrical structure. The stainless-steel armature goes through that. Another though hole perpendicular to that is tapped with M4 screw and the steel rod is held firmly using M4 setscrew. Material: Titanium Ti-6Al-4V
4. Wrld_Cam_Holder_Narrow.step	
	Description: This piece is the base for holding the world cam. The large hole through the short axis is where the stainless-steel armature passes through. It has a hole which I tap for m2 screw and use a metal m2 screw to hold the piece on to the armature. Material: Titanium Ti-6AI-4V
5. Wrld_Cam_Extension.step	
	Description: This piece holds the world camera. The 4 legs are glued to the 4 back corners of the camera sensor. The large through hole is meant for a long M2 screw to pass through and hold this piece to the "Wrld_Cam_Holder_Narrow" Material: Titanium Ti-6Al-4V

CEREBRO Head Assembly Parts Description Author: Dr. Vikram Pal Singh

Dept. of Psychology, UCSD Date: 01/06/23

Dept. of Psycho	logy, UCSD
6. Eye_Cam_Holder_Narrow.step	
	Description: This piece holds the eye camera. It has a through hole for the steel armature to pass through and a hole for M2 screw to hold the piece in place. The flat side of the piece is used to glue the eye cam while fitting on an animal. Material:
	Titanium Ti-6Al-4V
7. LED_Holder.step	
	Description: This piece is similar to other pieces with m2 screw hole to hold on to the steel armature and a large through hole for holding the IR LED that is used to illuminate the eye of the animal. It is always placed right next to the "Eye_Cam_Holder_Narrow"
	Material: Titanium Ti-6Al-4V
8. Mirror_Extension_Holder_V3.step	
	Description: This piece has a through hole to allow the "Mirror_Extension" to pass through which is held securely using a m4 setscrew. Material: Titanium Ti-6Al-4V
9. Mirror_Extension.step	
	Description: This piece is a key piece as it allows for adjustment of thee position of the hot mirror for the animal. The shaft goes through "Mirror_Extension_Holder_V3" and the placement is adjusted by moving the shaft up/down and can be rotated as well to create better position for the hot mirror Material: Titanium Ti-6Al-4V

CEREBRO Head Assembly Parts Description Author: Dr. Vikram Pal Singh

Dept. of Psychology, UCSD Date: 01/06/23

10. Hot_Mirror_Holder_Cube	
	Description: This cube holds the "Mirror_Extension" and "Hot_Mirror_Holder" together. All the holes are designed for M4 tap. Material: Titanium Ti-6Al-4V
11. Hot_Mirror_Holder.step	
	Description: This piece holds the hot mirror that is cut to size. Material: Titanium Ti-6Al-4V